

DISPOSABLE NONWOVEN BOOT COVER (ANTI-SLIP)

TECHNICAL DATA SHEET



CHARACTERISTIC	SPECIFICATION
Material Composition	Polypropylene Nonwoven Fabric, Non-slip rubber sole
Description	<p>Shoe Cover made of non-woven fabric</p> <p>32x38 cm (Suitable for all shoe sizes)</p> <p>Disposable Non Woven boot cover, lower protection device creating a real barrier against liquids and bacteria. Good tensile strength, both longitudinal and transverse.</p> <p>Made in TURKEY</p>
Basic Weight (g/m ²)	25 g/m ² (±2)
Colour	Blue /White